twinphoton EN/GR

User Manual

RED LIGHT THERAPY

User Manual

Create a better life with LED technology

Contents

Introduction.....	2
Red Light Therapy.....	3
Benefits.....	4
Parts Included.....	5
Instructions.....	6
Product Characteristics.....	7
How to use the Torch.....	8
FAQ.....	9
Storage & Safety.....	11

Introduction

Welcome to your new journey towards enhanced wellness and rejuvenation with the twinphoton Red Light Therapy Panel. Your decision to incorporate our advanced device into your daily wellness routine is a significant commitment to your personal wellness and vitality.

This revolutionary device uses high-quality LED lights to emit red and near-infrared light, which has been proven to provide many benefits for your personal well-being. The device's ease of use, safety features, and effective performance make it the ideal choice for your comprehensive wellness routine.

Before we start this wellness journey together, we kindly ask you to carefully read the usage instructions. This manual provides you information about the use, maintenance, and safety of your new device.

Red Light Therapy

What is Red Light Therapy?

Red light therapy, often referred to as Low-Level Light Therapy (LLLT) and Photobiomodulation (PBM), has multiple benefits, such as enhancing cellular function and supporting energy production system throughout the body. Twinphoton's red light therapy devices have numerous benefits for your body and have been proven to improve your sleep quality, skin appearance, relieve you from marks and scars, and assist in relaxing your muscles after exercise.

How does it work?

Certain frequencies of Red and Near Infra-Red light have beneficial properties and can enhance energy production in cells by stimulating proteins called cytochrome oxidase. This helps mitochondria break down nitric oxide and effectively produce more Adenosine Triphosphate (ATP). ATP is the cell's energy "fuel" and optimal ATP production supports and accelerates their restoration. It's important to note that most people in developed countries aren't exposed enough to the sun. Sun exposure is essential for human health and well-being. Studies have shown that 90% of our time is spent indoors. Red light therapy is a simple and quick solution to this problem, allowing cells to receive the necessary energy and care they need daily, regardless of weather conditions.

Benefits

Wrinkle Reduction: It can contribute to the reduction of the appearance of wrinkles, giving your skin a more youthful look.

Better Skin Quality: It can help reduce skin roughness, making it softer and more elastic.

Scar Appearance Reduction: It can improve the appearance of scars, making your skin more smooth.

Reduced Muscle Sensitivity: It can help reduce muscle sensitivity after exercise, speeding up recovery and return to normal activity.

Reduced Morning Stiffness: It can help reduce morning stiffness, improving your mobility and flexibility.

Better Sleep: It can contribute to improving the quality of sleep, enhancing the sense of rest and endurance ability throughout the day.

Anxiety Reduction: Red light therapy can help relieve anxiety, creating a sense of calm and relaxation.

Parts Included

x1 Red Light Therapy Torch

x1 Torch case

x2 Rechargeable batteries type 18650

x1 Battery Charger

x1 Protective Glasses

x1 User Manual

Operating Instructions


Wavelength:

- 630 nm
- 660 nm
- 850 nm

Timer: 5 minutes (Preselected): the device will automatically shut-off after 5 minutes

How to User the Torch

Examples of Using the Torch

Wrist

Fingers

Knee

Do not use the device continuously for more than 5 minutes to avoid overheating. If you want to use the device in more than one spots allow the device to cool for 1 minute between uses.

Product Characteristics

Product Size	126 x 30 mm	Weight	0.2 kg
Material	Stainless Steel	Wavelength	630: 660: 850nm
Irradiance	> 250mw/cm2@surface	Irradiance	8 cm distance 51mw/cm ²
Luminence	180 lux at 15cm	Battery	18650
Warranty	2 years	Lifespan	50.000 hours

Battery Placement

1. Place the included battery inside the charger
2. Place the positive pole in the **TOP** side(see picture)
3. Place the negative pole against the **BOTTOM**(see picture)
4. Use **only** 18650 batteries
5. Charge the **18650** batteries for 2-3 hours

Recommended Use

Red light therapy is used for wellness devices and is not intended for the diagnosis, prevention, treatment, or relief of disease, injury, or disability. The use of these devices helps promote healthy skin, muscle restoration, and wellbeing. If you have any questions or concerns about using these devices with a pre-diagnosed or underlying condition, please consult your doctor or healthcare provider before use.

Photosensitivity Test: It is important to perform a photosensitivity test before using the device for the first time. Activate the device in the red and near-infrared light emission mode and use it on your (bare) arm at a distance of 8 to 30 centimeters for 4 minutes. If you feel any discomfort or irritation on your skin, do not use the device unless you first consult a doctor or healthcare provider.

Protective Glasses: The emitted surface covers the entire front surface of the device. Glasses are included and should be used when using the device with the Near Infrared (NIR) Mode active. During treatment, do not look directly at the light-emitting diodes (LEDs).

First Use: If you are using the device for the first time, we recommend starting with 1-2 minute sessions per area, gradually increasing to 5-minute sessions over the course of 2-3 weeks as your body adjusts to the light exposure.

Session Time: As a rule, we suggest using it for 5 minutes per session, 15-30 centimeters away from your device, which is ideal. 5 minutes for symptomatic areas may be beneficial, but beyond that, the evidence shows that returns diminish (your cells only absorb so much light at a time). For this reason, we

recommend use once daily.

Frequently Asked Questions

How to use the Device

For maximizing cellular health, including a beneficial impact on deep tissues, it is recommended to use red light every day or at least five times a week, once daily.

Should I use protective glasses?

Protective glasses are necessary when your eyes are in direct contact with the light. Your eyes are not accustomed to red light, and the red light can seem very bright. For this reason, we recommend you to close your eyes or wear sunglasses. If you have any concerns, please consult your doctor.

Why do near-infrared LEDs seem like they are not working?

The human eye can only see light wavelengths up to an average of 760nm. Red light therapy uses wavelengths of 850nm for near-infrared LEDs, which are not visible to the naked human eye.

What is the difference between red and near-infrared light wavelengths? Both near-infrared light at 850nm and red light at 660nm are designed to enhance cellular function. The main difference is that skin tissues absorb red light, leading to increased collagen production and improved skin health. Near-infrared light at 850nm is not visible to the human eye. It penetrates lower levels in the tissues, benefits of which have been proven, such as muscle recovery.

Can I wear clothes during the use of the red light device?

Yes, you can wear clothes during the treatment. The only thing that matters is that the light should be able to reach your skin. The light cannot penetrate through clothes, so any area covered by clothing will not receive the benefits of the treatment. For maximum effectiveness, it's best to expose the treatment area directly to the light, without any clothing, makeup, or anything else that might block the light. However, it's important to protect sensitive areas of your body that may be more vulnerable to light exposure.

Will the red light interact with the medication I'm taking?

Many studies have shown that phototherapy is safe and effective. However, if you are using photosensitizing drugs, steroid medications, or injections, you should consult your doctor before using the device. If you feel any discomfort or concern while using the device, contact your doctor immediately.

Is it safe for children and pregnant women to use the device?

There are no data for red light therapy pregnant women and children. Therefore, its use is not recommended.

Storage and Safety

Before cleaning your device(s), make sure the device(s) is turned off and the power cords are disconnected from any power source. Before cleaning, remove the devices from their mounting accessories. To clean your panel, except for the lens parts, you can use a microfiber cloth with water to clean the body of the device. Do not spray or apply moisture directly to the device. Be careful when wiping around the openings of ventilation slots, LED openings, and any input or output openings to avoid moisture ingress. ***Caution: Do not use chemical cleaners, such as alcohol, to clean the lens surface of the device.**

Storing the Device

The device and accessories should be stored in a clean and dry place.

***Temperature 0°C ~ 40°C (32°F -104°F)**

***Humidity less than 50%**

If the device presents any problem, please contact us. The devices and the device's components cannot be serviced by the user.

The devices and components are covered by the two-year statutory warranty of Article 5 of Law 2251/1994 and 534 of the Greek Civil Code.

Warning

1. The device should not be used by individuals with reduced physical, sensory, or mental capabilities.
2. Do not allow children to play with the device.
3. **WARNING:** Do not stare into the emitter with your eyes. Read the instructions carefully.
4. Please protect your eyes from exposure to visible and infrared radiation. Strictly follow the cycle time provided in the usage instructions
5. If unexpected side effects, such as itching, occur within 48 hours of first using the device, medical advice should be sought before further exposure.
6. The device should not be used if the timer is faulty or if the filter is broken or has been removed. If you notice any damage to the device, do not use it.
7. User replacement is prohibited, special tools are required for the replacement of the LED diodes.
8. Do not use the device near water or liquids.
9. Do not put the device in running water.
10. **If the device is not going to be used for more than 1 month, remove the batteries from the device to avoid possible damage.**

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